



**EXPEDITED SERVICE
FOR LIQUID HANDLING DEVICES ONLY**
 YES NO
 See Page 2 for Details

Ship To: Calibration Lab
 BioTek Services, Inc
 5310 S Laburnum Ave
 Henrico, VA 23231
 (804) 222-5833
www.biotekcalibrations.com

Customer Shipping Info:

Company:		Contact:	
Street		Dept:	
City, State, Zip:		Rm/lab:	
Email:		Phone:	

Billing Info: Check here if billing is the same as shipping

Company:			
Street			
City, State, Zip:		Phone:	
Contact:		Fax:	

Payment Options:

Purchase Order #			
Credit Card #:			
Expiration Date:		CVC Code	
Signature approving charge			

Return Shipping: (Check One) Overnight 2 days Ground
 UPS _____ FedEx (account number) _____

Pipettes:	Price	Quantity		Price	Quantity
Level 1 Single Channel	\$33		Level 1 Multichannel (8 & 12) Channel	\$125	
Level 2 Single Channel	\$28		Level 2 Multichannel (8 & 12) Channel	\$95	
Level 3 Single Channel	\$18		Level 3 Multichannel (8 & 12) Channel	\$45	

Calibration frequency: 3 months 6 months 12 months Other: _____

BottleToppers/Dispensettes:

	Price	Quantity		Price	Quantity		Price	Quantity
Level 1	\$75		Level 2	\$50		Level 3	\$25	

Calibration frequency: 3 months 6 months 12 months Other: _____

All Instrument Must Be Decontaminated OR Sterilized Prior To Shipping

Autoclave Warning: Autoclaving may damage your equipment, please verify with manufacturer specifications for autoclavability.

Standard test points are typically: 100%, 50%, and 10% of the nominal with reads at high, medium, and low unless otherwise specified on the pipette.

NOTE: Please include several tips commonly used in your lab with your pipettes, for validation purposes.

Level One service includes (ISO 17025): PM, Calibration and "AS-Found" and "As-Left" Data Reports

Level Two service includes (ISO 17025): PM, Calibration and "As-Left" Data Reports

Level Three service includes (NIST Traceable): PM, Calibration and a Certificate of Calibration.

Temperature:			Quantity	Quantity		
Liquid in glass Thermometers Must be shipped ground (ISO) Each additional Test Point \$20	1 point	\$50		3 pts	\$90	
Temperature Device with Probe Each additional probe (ISO)	1 Point	\$75 \$25		3 pts	\$115 \$65	
Thermo-Hygro Temp & Humidity (15% - 95%) NIST Traceable	1 Point	\$80		3 pts	\$ 120	
Thermo-Hydro Temp, Humidity & Bar Pressure NIST Traceable	1 Point	\$85		3 pts	\$ 125	

Calibration frequency: 3 months 6 months 12 months Other: _____

Other Calibrations:	Pricing	Quantity	w/battery change \$30
Timers & Stopwatches	\$20		
pH-Measuring Equip. (4,7,10)	\$45		

Calibration frequency: 3 months 6 months 12 months Other: _____

Special Notes or Instructions

Disclaimer:

This is not an absolute guarantee and BioTek Services, Inc. will use all reasonable efforts to meet the expressed time for performance. BioTek Services, Inc. reserves the right to void this service offer or adjust the fees and other charges for any reason including large quantities of instruments. Customer will be notified of all exceptions, changes, or increase in service fees using every reasonable means of communication. The Company's liability to the Client or any third party, for a claim of any kind arising as a result of or related to any product or service provided by BioTek Services, Inc., whether in contract, in tort (including negligence or strict liability) or otherwise, under any warranty, condition or guarantee or otherwise, shall be limited to the cost of this service or \$75.00, whichever is greater. BioTek Services, Inc. assumes no liability or responsibility for damages or errors incurred during shipping or transit. Such damages or errors must be claimed against the carrier. Due to a high rate of lost, misdirected, and delayed deliveries, utilizing the U.S. Postal Service for package delivery is not recommended.

We follow simple acceptance, statements of conformity (e.g. Pass/Fail) to specifications are made in this report without taking measurement uncertainties into account. However, measurement uncertainties are provided for your utilization. Pass= Results within limits/specifications Fail= Results exceed limits/specifications

Turnaround Time:

The turnaround time is three-five business days. "Turnaround time" is defined as the actual or realized chronological length of service time in our facility excluding transportation or shipping time.

EXPEDITED SERVICE: (For Liquid Handling Devices Only)

We can provide **EXPEDITED SERVICE** on most calibrations for an additional charge of \$75.00 for each level.

Fee Limits for Single Channel Pipettes: 50 Level Three, or 25 Level Twos, or 15 Level Ones.

Fee Limits for Multi-Channel Pipettes: 4 (12 Channel Pipettes) or 6 (8 Channel Pipettes), Level I

Fee Limits for Multi-Channel Pipettes: 6 (12 Channel Pipettes) or 8 (8 Channel Pipettes), Level II

Fee Limits for Multi-Channel Pipettes: 8 (12 Channel Pipettes) or 10 (8 Channel Pipettes), Level III

Fee Limits for Bottle Toppers/Dispensettes: 5 Level I, or 8 Level II, or 10 Level III

The turnaround time for EXPEDITED SERVICE is one business day, excluding weekends and legal holidays, provided the items are received in our laboratory before 2:00pm (EST). The shipping fees are not included in the Expedited Service Fee.

Services Performed Onsite Only:** Scales, Balances, Microscopes, Centrifuges, Ovens, Refrigerators, Freezers, Coolers, water baths & Incubators. ***Call for Pricing